



PRESTAN®

AED UltraTrainer

Instructions for using the Programming Dongle



Basic steps:

1. Make sure the unit has good batteries in it, and unit is off
2. Turn unit over and remove gray rubber plug from socket as shown in the pictures at right
3. Insert programming dongle into the socket, aligning beveled corners on dongle with beveled corners in socket – push down until dongle is flush with unit
4. Turn unit back over and set on table top with dongle tab dangling over the table edge
5. Press green power button (this will turn on the unit and run the programming sequence for approximately 2 minutes / lights will flash during the programming)
6. When programming is done, unit will beep and shut off
7. Press green power to turn on unit, you should hear “Adult Patient...”
8. (If you don’t hear “Adult Patient...” program again)
9. Press green power button and hold down for 3 seconds to turn unit off
10. Remove programming dongle
11. Insert gray plug back into socket
12. Unit is ready to use

