

Practi-TRAINER Essentials by



Available at www.MCRmedical.com



Practi-TRAINER™ Essentials is affordable and compact enough for instructors to bring one for everybody in the class! It offers a universal training experience to prepare students to use a variety of different AEDs. It features two preprogrammed scenarios, clear voice prompts, built-in Bi-Lingual options, and a metronome. **Practi-TRAINER™ Essentials** is compatible with any type of CPR manikin that is intended for AED training.

Key features and benefits are realistic defibrillator training, and affordable enough for each student in a class to have their own trainer. Voice prompts guide students through simulated rescue situations. Metronome sets proper pace for effective chest compressions. Lightweight, compact unit is ideal for transporting to trainings at any location.

Each unit comes with:

- AED trainer
- 1 set of reusable adult training pads and connector
- 1 set of reusable child training pads and connector
- 3 AA batteries
- Nylon carrying case
- Instructor manual
- 1-year manufacturer's warranty

Individual trainer measures 5" x 5" x 1.75"

Weight: 6.9 oz (trainer only)

Single unit item# **WL120ES10**

Case measures 9" x 7" x 2.25"

4-Pack item# **WL120ES10-4**

Case measures 14.25" x 14.25" x 3.5"



MCR | **MEDICAL**

\$69.95 (WL120ES10) / \$249.95 (WL120ES10-4) at
www.mcrmedical.com/category/AED-Trainers.html