



**PRESTAN<sup>®</sup>**

# PRODUCT CATALOG

Available at  
**MCRmedical.com**



# PRESTAN® PROFESSIONAL AED TRAINER

## Prestan® Professional AED Trainer

The **Prestan Professional AED Trainer** provides quality, durability and affordability like no other trainer in the market. Clear, confident voice instructions that provide easy-to-follow directions and preconfigured scenarios make training students easier. The Prestan Professional AED Trainer carries a three (3) year warranty.



- Easy to use panel includes a scenario, volume control, and language button
- Includes 5 scenarios for varied training opportunities
- Pads can be used at least 50-60 times during CPR/AED classes
- Pads are pre-connected
- Pads have an embedded 'switch' which senses proper pad attachment to the manikin
- The silicone base adhesive on the pads, coupled with the foam structure, ensures easy application and removal with no residue left behind
- Pads carry a one (1) year warranty
- Module contains 2 languages (English / Spanish)\*
- Module contains the most up to date guidelines and language capability
- When guidelines change, only the module needs to be replaced



### Optional Remote Control

- Volume, language and scenario control
- Simulate loose pad
- Initiate "Press Deeper" during CPR
- Pause, play or stop



### Optional Pediatric Pads



Available at  
**MCRmedical.com**



# PRESTAN® PROFESSIONAL AED TRAINER



PP-AEDT-101

## Prestan Professional Single AED Trainer

- Includes one set of Adult Pads w/Case, batteries, instructions and a carrying case
- Compatible with any CPR manikins
- Can be used to train in both semi-automatic and fully-automatic versions
- Metronome and counts upward by 10's ending at 30 to help the student with chest compressions
- KIT includes Pediatric Pad w/Case and Remote Control

## **PRESTAN AED TRAINER ITEM LIST**

- # PP-AEDT-101 (Single AED Trainer)
- # PP-AEDT-401 (4-Pack AED Trainers)
- # PP-APAD-1 (Single Set of Adult Pads)
- # PP-APAD-4 (4-Pack Sets of Adult Pads)
- # PP-ACASE-1 (Single Adult Pads Case)
  
- # PP-AEDT-KIT-101 (Single AED Trainer Kit)
- # PP-AEDT-KIT-401 (4-Pack AED Trainer Kit)
  
- # PP-PPAD-1 (One Set of Pediatric Pads)
- # PP-PPAD-4 (Four Sets of Pediatric Pads)
- # PP-PCASE-1 (One Pediatric Pads Case)
- # PP-PPDC-1 (One Set of Pediatric Pads w/Case)
  
- #PP-AEDT-100-R (Single Remote)
- #PP-AEDT-400-R (4-Pack Remotes)
- ITEM #PP-AEDT-101-M (Single Module)
- ITEM #PP-AEDT-401-M (4-Pack Modules)

*\* Other Languages Available*



PP-AEDT-KIT-101



PP-AEDT-KIT-401



PP-AEDT-401



Available at  
**MCRmedical.com**



# PRESTAN® PROFESSIONAL TRAINING MANIKINS

## Prestan® Professional CPR-AED Training Manikins

Realistic to the eye and the touch, the **Prestan Professional Training Manikin** series is unlike any other on the market. Available individually or in convenient multi-packs of four for class training purposes, our manikins are uniquely designed to give each student a realistic experience that provides for the best "real life" CPR outcomes. Prestan Manikins carry a three (3) year warranty.

### MANIKIN SINGLES



Manikin single units come with  
10 Face-Shield/Lung-Bags  
and a nylon Carrying Case

- Available in Adult, Child and Infant sizes
- Light (original), Medium and Dark skin tones available
- Realistic looking and performing manikin is more acceptable to both instructors and students.
- Far easier and faster to carry and setup at training sites
- Easily inserted face-shield/lung-bag provides visual chest rise to ensure students receive the most accurate and realistic training
- Face shield is secured so it won't slip and there's no need for adhesives
- Patented face/head tilt simulates how an actual victim's head would move if he required CPR and teaches students how to properly open the airway
- Manikins are available with or without the CPR Feedback Monitor - Monitors may also be purchased separately to easily install into non-monitored Prestan Professional Manikins
- Prestan Manikins are latex free

### MANIKIN 4-PACKS



Shown in Light  
Skin Tone



Shown in Medium  
Skin Tone



Shown in Dark  
Skin Tone

Manikin 4-packs come with  
50 Face-Shield/Lung-Bags  
and a nylon Carrying Case



Available at  
**MCRmedical.com**





# PRESTAN® PROFESSIONAL TRAINING MANIKINS

## AMAZING FEEDBACK TOOLS MAKE ALL THE DIFFERENCE!

Prestan's revolutionary *CPR Rate Monitor* allows for instant feedback to both instructor and student regarding the rate of chest compression. This allows each student to gauge their rate of compressions on their own as well as allowing the instructor to monitor several students quickly and easily.



The unique *Clicker Mechanism* provides an audible confirmation when the correct compression depth is reached. This increases student confidence by guiding to the appropriate force when administering CPR compressions.

## ADULT JAW THRUST MANIKIN

Prestan's Jaw Thrust Head gives instructors the option to teach with the jaw thrust maneuver or head-tilt/chin-lift.



Prestan Professional  
**Adult Jaw Thrust**  
CPR-AED Training Manikin

*Jaw Thrust Heads may also be purchased separately - they are interchangeable with the heads on the Prestan Professional Adult Training Manikins*

## MANIKIN COLLECTIONS

Prestan Manikins are available in 2 convenient multi-packs.



**The Prestan Collection**  
(1 Adult / 1 Child / 1 Infant)  
Comes with 10 Face-Shield/Lung-Bags  
per Manikin



**The Prestan Family Pack**  
(2 Adults / 1 Child / 2 Infants)  
Comes with 10 Face-Shield/Lung-Bags  
per Manikin




Available at  
**MCRmedical.com**



# PRESTAN® ULTRALITE® MANIKIN 4-PACK

## Prestan® Ultralite® Manikins

Prestan is pleased to offer our newest and most portable manikin, the **Ultralite Manikin**. The Ultralite Manikin is amazingly easy to set up and is sold in a convenient and lightweight 4-pack for efficient training-on-the-go. This durable manikin offers an affordable tool for CPR training, with all the quality and realism you expect from Prestan.



- Uniquely fast and easy set-up
- Stackable torsos and heads generate a compact package for transporting multiple units
- Clamshell opening provides easy access
- Manikin is easy to clean
- Vinyl and latex free
- Compatible with AED Trainer pads
- Realistic head tilt
- Visible chest rise with lung bag use
- Unique torso web design for realistic feel
- Anatomically correct design includes xyphoid process
- Realistic chest compressions to a 2" depth
- Life tested to a half-million compressions
- Three (3) year warranty

### Ultralite Manikin 4-pack includes:

- 4 Manikin Torsos
- 4 Manikin Heads
- 4 Compression Pistons
- Nylon Piston Bag
- 50-pack Face-Shield Lung Bags
- Nylon Carrying Case
- Detachable Shoulder Strap
- Instruction Sheet



Prestan Ultralite Manikins  
ITEM # PP-ULM-400-MS



**Complete  
4-pack  
weighs only  
13 pounds!**



Available at  
**MCRmedical.com**



# PRESTAN® PRODUCTS - ACCESSORIES

## Prestan® Products - Accessories

### DON'T RUN OUT!

Plan ahead to be sure you have enough Prestan accessories for your training needs.  
Prestan Accessories are available from your Prestan Authorized Distributor.



**Prestan Professional  
Face-Shield Lung Bags**  
Adult 50-Pack: PP-ALB-50  
Adult 10-Pack: PP-ALB-10  
Also available for  
Child, Infant and Ultralite Manikins



**Prestan Professional  
Face Shields**  
Adult 50-Pack: PP-AFS-50  
Adult 10-Pack: PP-AFS-10  
Also available for  
Child, Infant and Ultralite Manikins



**Prestan Professional Rescue  
Mask Adaptors**  
(10-Bag)  
10076-PPA



**Prestan Professional  
Generic Face Shield**  
Adult 36-Count  
10077-PPS

*Accessories are considered consumable items meant for a single use and are not covered by  
Prestan's three (3) year product warranty.*



Available at  
**MCRmedical.com**





Prestan Products, LLC  
701 Beta Drive, Suite 3 • Mayfield Village, OH 44143 • 800-434-7989  
[www.prestanproducts.com](http://www.prestanproducts.com)

© January, 2015

Available at  
**MCRmedical.com**

