



Practi-CRdM® Specification Sheet

Item Number:	WLCCRdM	Case Pack:	20
Size	Case: 3."x2.75"x2 – Band 10"x0.75"	Weight:	0.1125 lbs
Series:	CPR Training	Case Weight:	2.25 lbs
Customization:	N/A	Case Dimensions:	17.5"x7"x5"

Description	Image
<p>Practi-CRdM Wristband</p> <p><i>The NEW Practi-CRdM® Compression Rate and depth Monitor is designed to measure compression rate and depth when practicing CPR on a manikin. Worn on the students' wrist, a diode illuminates red when compressions are below 100 cpm or above 120 cpm. When compression rate is between 100 and 120 cpm the diode illuminates green, allowing users to self-check their timing. In addition, the CRdM will provide an audible "beep" when compression depth is not within the range of 5-6cm (2"-2.3").</i></p> <p><i>A flexible wristband fits any size wrist and may be used with any training manikin. LATEX FREE</i></p> <p><i>The Practi-CRdM® provides an affordable solution that helps protect the investment training centers have in the manikins they own which do not offer compression rate and or depth monitoring, key elements the AHA is requiring by January 31, 2019.</i></p>	 <p style="font-size: 2em; font-weight: bold; margin-top: 10px;">Practi-CRdM®</p>

Features & Content

Features:	<ul style="list-style-type: none"> Meets the January 2019 AHA mandate regarding the use of an "instrumented directive feedback device" Easy to use; little or no special training of instructors needed Works with any manikin! Measures BOTH compression rate and compression depth Adjustable wristband fits nearly any size wrist Economical Fraction of the cost of buying new manikins or retrofitting existing manikins Rechargeable lithium battery voltage of 3.7V and capacity of 90mAh via USB port Minimum 3-4 hours of <u>continuous</u> use on one battery charge One year warranty 	Packaging:	<p>■ PANTONE Black C</p> 
-----------	---	------------	--

More Information

Complies With The Following Regulations:	CE
Manufactured For:	WNL Products
Purchasing Manager:	JF
Assembly:	Assembled in China

Last Updated: 1.22.2019

