

# Which Prestan® AED Trainer is Right for Me?



## >AEDT Plus vs. UltraTrainer™<

In the market for a new AED Trainer? Which one best suits your training needs?



<u>Prestan Professional AEDT Plus</u>	<u>Features</u>	<u>Prestan AED UltraTrainer</u>
Similar to Full-Size AED Unit is 10" x 7.25" x 2.75" and weighs 2.35 lbs.	<b>Size and Weight</b>	Compact and Portable Unit is 7.25" x 5.25" x 2" and weighs 1.2 lbs.
Yes	<b>CPR Feedback</b>	Yes
Yes	<b>Pad Sensing System (PSS)</b>	Yes
No	<b>Pre-connected Pads</b>	Yes
5 Scenarios Available	<b>Shock Simulation</b>	5 Scenarios Available
Yes	<b>Auto / Semi-Automatic Shock</b>	Yes
Yes	<b>Rescue Breaths</b>	Toggle On/Off
Yes	<b>Compression Metronome</b>	Toggle On/Off
Yes	<b>Child Button</b>	Yes
Dual Adult/Child Pads	<b>Pediatric Pads</b>	Dual Adult/Child Pads
Available	<b>Remote Control</b>	No
Included	<b>Long-Life "C" Batteries</b>	Included
Yes	<b>Low Battery Prompt</b>	Yes
3 Years	<b>Prestan Warranty</b>	3 Years

Visit [www.prestanproducts.com](http://www.prestanproducts.com) for full details on these quality products.  
Contact your Prestan Authorized Distributor for pricing and ordering.

