


Prestan® Ultralite™ Manikins

Prestan is pleased to offer our newest and most portable manikin, the **Ultralite Manikin**. The Ultralite Manikin is amazingly easy to set up and is sold in a convenient and lightweight 4-pack for efficient training-on-the-go. This durable manikin offers an affordable tool for CPR training, with all the quality and realism you expect from Prestan.



Ultralite Manikins 4-pack features:

- Uniquely fast and easy set-up
- Stackable torsos and heads generate a compact package for transporting multiple units
- Clamshell opening provides easy access
- Manikin is easy to clean
- Vinyl and latex free
- Compatible with AED Trainer pads
- Anatomically correct design includes xiphoid process
- Realistic head tilt
- Visible chest rise with lung bag use
- Unique torso web design for realistic feel
- Realistic chest compressions to a 2" depth
- Life tested to a half-million compressions
- Three (3) year warranty

Ultralite Manikins 4-pack includes:

- 4 Manikin Torsos
- 4 Manikin Heads
- 4 Compression Pistons
- Nylon Piston Bag
- 50-pack Face-Shield Lung Bags
- Nylon Carrying Case
- Detachable Shoulder Strap
- Instruction Sheet

Reorder from your distributor:

- **#PP-ULB-50**
(50-pack Ultralite Manikin face-shield/lung-bags)
- **#PP-UFS-50**
(50-pack Ultralite Manikin face-shields)



ITEM # PP-ULM-400-MS
(Prestan Ultralite Manikins)



Complete 4-pack weighs only 13 pounds!